Molecule Set Comparator (MSC)

- A short introduction -

Jan-Mathis Hein






Starting the application

- Open the "MSC_1.0" folder.
- Select and execute one of "MSC_?GB_<OS>" files where <OS> should be your operating system and? specifies the amount of RAM that must be allocated for the JVM in gigabytes.

lame ^	Änderungsdatum	Тур	Größ
GPLv3.txt	31.05.2020 09:11	Textdokument	
MSC_0.7_Tutorial.pdf	31.05.2020 09:11	Adobe Acrobat D	
MSC_0.7_Tutorial.pptx	31.05.2020 09:11	Microsoft PowerPo	5
MSC1GB_Linux.sh	31.05.2020 09:11	SH-Datei	
MSC1GB_Mac.sh	31.05.2020 09:11	SH-Datei	
MSC1GB_Windows.bat	31.05.2020 09:11	Windows-Batchda	
MSC2GB_Linux.sh	31.05.2020 09:11	SH-Datei	
MSC2GB_Mac.sh	31.05.2020 09:11	SH-Datei	
MSC2GB_Windows.bat	31.05.2020 09:11	Windows-Batchda	
MSC4GB_Linux.sh	31.05.2020 09:11	SH-Datei	
MSC4GB_Mac.sh	31.05.2020 09:11	SH-Datei	
MSC4GB_Windows.bat	31.05.2020 09:11	Windows-Batchda	
MSC8GB_Linux.sh	31.05.2020 09:11	SH-Datei	
MSC8GB_Mac.sh	31.05.2020 09:11	SH-Datei	
MSC8GB_Windows.bat	31.05.2020 09:11	Windows-Batchda	
MSC16GB_Linux.sh	31.05.2020 09:11	SH-Datei	
MSC16GB_Mac.sh	31.05.2020 09:11	SH-Datei	
MSC16GB_Windows.bat	31.05.2020 09:11	Windows-Batchda	
MSC32GB_Linux.sh	31.05.2020 09:11	SH-Datei	
MSC32GB_Mac.sh	31.05.2020 09:11	SH-Datei	
MSC32GB_Windows.bat	31.05.2020 09:11	Windows-Batchda	
VersionHistory.txt	31.05.2020 09:11	Textdokument	

Initializing a new job


- Create a new job by choosing Job
 New in the menu bar.
- Then click on **Choose File** button.
- Navigate to the "TestSets" folder under the MSC folder and select the "testSet1.txt" file.
- Repeat the same steps to select the "testSet2.txt" file.
- Finally, press the Tanimoto and Atom counts drop down menu (marked with red arrows) to expand them.






Initializing a new job

- Select the Basic fingerprint and Atom count checkboxes.
- Then start the job by pressing the **Start** button and wait until the job gets over.



The output tab

- You can save the charts as PDF files by clicking on the Image button.
- Switching to the next descriptor can be done using the drop down menu.
- To get more fine detailed information you can click on the bins.



Center

Default

Rel. Frea.

Image

55,000

50,000

45.000 -

40,000 -

35,000

30,000 -

25,000

20,000 -

15,000 -

5.000 -

151 117 84 87

Detail window

- Use the arrow buttons to browse through the molecule pairs or use the text field to go to a specific position in the list.
- To save a list of molecules that contains the molecule in the visualized in the left side, click on the L list button.
- For more information, additionally you can select the PubChem fingerprint checkbox to calculate another Tanimoto coefficient.
- To get back to the previous window, close the "Molecule pairs of bin 1" window.






Chart configuration


- Change the upper border of the last bin by using the marked slider
- Change the upper bound of the y-axis range by using the marked slider
- The bin borders and y-axis bounds can be set explicitly by pressing the **Configure** button can change the number of bins to 5 by typing in the text field and hitting the enter key
- The bin borders can also be set explicitly by pressing the Bin borders button





Saving the output

- To save the output ,select Job ->
 Save in the menu bar.
- Select or create a folder where you need the output files to be saved.
- These results can be loaded back by selecting Job -> Load in the menu bar and selecting the corresponding folder.



Changing the preferences

- Open the "Other preferences" dialog by choosing Preferences ->
 Others in the menu bar.
- Here,
 - The preferences for the number of parallel threads to be used for computations
 - The default number of bins.
 - The number of molecule pairs to be saved when we save a job .
 - And the image quality can be set.

