



Molecule Set Comparator (MSC)

- A short introduction -

Jan-Mathis Hein













**Westfälische
Hochschule**



**FRIEDRICH-SCHILLER-
UNIVERSITÄT
JENA**




Starting the application

- Open the “MSC_1.0” folder
- Select and execute one of the “MSC_<OS>” files where <OS> should be your operating system.
- To have a smooth performance, the RAM allocation to the JVM is set to 4GB as a default.
- If you want to start the MSC with a different RAM specification, you must edit one of the existing “MSC_<OS>” files by changing the “-Xmx?g” option in the third row. ? specifies the amount of RAM (in gigabytes) that will be allocated to the JVM

 jdk-11.0.2	11.07.2020 20:40	Dateiordner	
 lib	11.07.2020 20:40	Dateiordner	
 Licenses	11.07.2020 20:40	Dateiordner	
 TestSets	11.07.2020 20:40	Dateiordner	
 Tutorials	11.07.2020 20:40	Dateiordner	
 MSC_Linux.sh	11.07.2020 20:40	SH-Datei	1 KB
 MSC_Mac.sh	11.07.2020 20:40	SH-Datei	1 KB
 MSC_Windows.bat	11.07.2020 20:40	Windows-Batchda...	1 KB
 README.MD	11.07.2020 20:40	MD-Datei	2 KB
 VersionHistory.txt	11.07.2020 20:40	Textdokument	1 KB


Initializing a new job

- Create a new job by choosing **Job** -> **New** in the menu bar.
- Then click on the **Choose File** button.
- Navigate to the “*TestSets*” folder in the MSC folder and select the “*testSet1.txt*” file.
- Repeat the same steps to select the “*testSet2.txt*” file.
- Finally, click the **Tanimoto** and **Atom counts** drop-down menus (marked with red arrows) to expand them.




Initializing a new job

- Select the **Basic fingerprint** and **Atom count** checkboxes.
- Then start the job by pressing the **Start** button and wait until the job finished.



The output tab

- Once the execution completes, a separate tab with the results will pop up.
- You can save the chart as a PDF file by clicking on the **Image** button.
- Switch to the chart of the “Atom count” descriptor by pressing the upper choice box and choosing “Atom count”



Detail window

- To open the detail window, click on the marked bin.
- Use the arrow buttons to browse through the molecule pairs contained in this bin or use the text field to go to a specific position in the list.
- Save a list that contains the left molecule of every pair by pressing the **L list** button
- Select the **PubChem** fingerprint checkbox to calculate an additional Tanimoto coefficient
- Then close the “*Molecule pairs of bin 1*” window.






Chart configuration

- Change the upper border of the last bin by using the provided slider.
- Change the range of the y-axis by using the provided slider.
- The bin borders and y-axis range can also be set explicitly by pressing the **Configure** or **Bin borders** buttons
- Change the number of bins to 5 by typing in the text field and hitting the enter key



Saving the output

- To save the output, select **Job** -> **Save** in the menu bar.
- Select or create a folder where you want the output files to be saved.
- These results can be reloaded by selecting **Job** -> **Load** in the menu bar and selecting the corresponding folder.



Changing the preferences

- Open the “*Other preferences*” dialog by choosing **Preferences** -> **Others** in the menu bar.
- Here:
 - The preferences for the number of parallel threads to be used for computations,
 - The default number of bins,
 - The number of molecule pairs to be saved when saving a job,
 - And the image quality for saved images can be set.

